



ACCEPTING ALL DIGITAL SCAN FILES



To send a file:

1. Create an account on cerec-connect.com and follow the prompts to activate an account
2. Log into your account and select Pan-Am Dental Lab from the list of available labs
3. Add Pan-Am Dental Lab to your list of "my favorite Laboratories"
4. Complete your scan and online Rx form. Choose Model-Less if you do not want a model.
5. Select Pan-Am Dental to be the file recipient.



3M™ True Definition
Scanner

To send a file:

1. Log into 3M connection center. If you need assistance to log in call 877-722-6528.
2. Visit the Scanner Laboratory Menu and select Pan-Am Dental Lab.
3. Complete the scan and online Rx form. Make sure to indicate if you do not want a model printed.
4. Upload scan data to the 3M communication center.



iTero™

To Send a File:

1. Sign into the Cadent iTero account online and add Pan-Am Dental Lab to your scanner menu. If you have trouble contact 800-577-8767
2. Complete your scan and fill out the online Rx form. Be specific in the notes section if you do not want a model or other instructions.
3. Select Pan-Am Dental Lab to receive the file.



TRIOS®

To Send a File:

1. Go Online and register with 3Shape Communicate and send Connection Request to Pan-Am Dental Lab.
2. Pan-Am will respond that it accepts the request to send files.
3. Complete the intraoral scan and Trios Rx form. In the notes section be specific about what materials you want to use and any other information.
4. Select Pan-Am Dental Lab and Send File