

Implant	Torque (Ncm)	Torque Settings	Driver & Head
ASTRA TECH			
Osseospeed™ 3.5, 4.0	20 Ncm	0.050" Hex	⬡
Osseospeed™ 4.5, 5.0	20 Ncm	0.050" Hex	⬡
BIOHORIZONS			
Internal 3.5	30 Ncm	0.050" (1.25 mm) Hex	⬡
Internal 4.5	30 Ncm	0.050" (1.25 mm) Hex	⬡
Internal 5.7	30 Ncm	0.050" (1.25 mm) Hex	⬡
BIOMET			
Certain™ 3.4	20 Ncm	0.048" Hex	⬡
Certain™ 4.1, 5.0, 6.0	20 Ncm	0.048" Hex	⬡
HIOSEN			
Conic TSIII mini	20 Ncm	.050 (1.25 mm) Hex	⬡
Conic Regular	30 Ncm	.050 (1.25 mm) Hex	⬡
NOBEL BIOCARE®			
NobelActive™ NP	15 Ncm	Nobel Biocare® Unigrip	★
NobelActive™ RP	35 Ncm	Nobel Biocare® Unigrip	★
NobelReplace™ 6.0	35 Ncm	Nobel Biocare® Unigrip	★
NobelReplace™ NP	35 Ncm	Nobel Biocare® Unigrip	★
NobelReplace™ RP	35 Ncm	Nobel Biocare® Unigrip	★
NobelReplace™ WP	35 Ncm	Nobel Biocare® Unigrip	★
STRAUMANN®			
Bone Level NC (Narrow CrossFit™)	35 Ncm	SCS Screwdriver	★
Bone Level RC (Regular CrossFit™)	35 Ncm	SCS Screwdriver	★
ZIMMER DENTAL			
Tapered Screw-Vent 3.5	30 Ncm	1.25 mm Hex	⬡
Tapered Screw-Vent 4.5	30 Ncm	1.25 mm Hex	⬡
Tapered Screw-Vent 5.7	30 Ncm	1.25 mm Hex	⬡